



PAGE 1 OF 1 - REPORT 1807223

REPORT DATE: July 31, 2018

REPORT TO: Eva Lau Mayway Corporation 1338 Mandela Parkway Oakland, CA 94607 Tel: 510-208-3113 Email: evalau@mayway.comSUBJECT: PO Number: LZ xiao yao; Samples for Gluten Analysis; RUSH ORDER RECEIVED: July 30, 2018; UPS Next Day Air TRK# 12 89E 637 01 4761 1537

ANALYTICAL RESULTS:

CONTROL #	SAMPLE DESCRIPTION	ANALYSIS IDENTIFICATION	Gluten
1807223-1	MW#3675 Free & Easy Wanderer Teapills/ Xiao Yao Wan Batch# A7L121	500821 Quantitative Gluten ELISA	32.6 ppm

RESULTS APPLY ONLY TO THE SAMPLES TESTED

NOTES:

The ALLER-TEK Gluten assay utilizes ELISA techniques and is designed for the measurement of cereal gluten at low levels, nominally within a range of 5-80 PPM, in raw materials and finished food products; it detects the gliadin and glutenin fractions of gluten.

For the quantitation of barley gluten, a standard curve derived from barley proteins is used.

The ALLER-TEK Gluten ELISA has been validated and certified with license no. 081202 as a Performance Tested Method by the AOAC Research Institute.

Respectfully Submitted for ELISA Technologies, Inc.

By _____ Analyst

Ву

Emily Warnock, Testing Supervisor

End of Report

THIS REPORT SHALL NOT BE REPRODUCED EXCEPT IN ITS ENTIRETY AND WITH THE WRITTEN PERMISSION OF ELISA TECHNOLOGIES, INC.

2501 NW 66th Court, Gainesville, Florida 32653, USA Telephone: (352) 337-3929 Fax: (352) 337-3928 E-mail: info@elisa-tek.com