



ELISA Technologies, Inc.

www.elisa-tek.com



PAGE 1 OF 1 - REPORT 1809212

REPORT DATE: October 2, 2018

REPORT TO: Eva Lau
Mayway Corp.
1338 Mandela Parkway
Oakland, CA 94607
Tel: 510-208-3113 ext. 8131
Email: evalau@mayway.com

SUBJECT: PO Number: QC; Samples for Gluten Analysis

RECEIVED: September 27, 2018; UPS Ground TRK# 1Z 89E 637 03 4645 9399

ANALYTICAL RESULTS:

CONTROL #	SAMPLE DESCRIPTION	ANALYSIS IDENTIFICATION	Gluten
1809212-1	Great Invigorator Teapills/Shen Tong Zhu Yu Wan Batch#2016012 MW#3390	500821 Quantitative Gluten ELISA	Not Detected < 2.5 ppm
████████	██ ██	████████████████████ ████████████████	████████████████ ████████████████

RESULTS APPLY ONLY TO THE SAMPLES TESTED

NOTES:

The ALLER-TEK Gluten assay utilizes ELISA techniques and is designed for the measurement of cereal gluten at low levels, nominally within a range of 2.5-80 PPM, in raw materials and finished food products; it detects the gliadin and glutenin fractions of gluten.

The ALLER-TEK Gluten ELISA has been validated and certified with license no. 081202 as a Performance Tested Method by the AOAC Research Institute.

Respectfully Submitted for **ELISA Technologies, Inc.**

By Tom d'Angelo
Analyst

By Emily Warnock
Emily Warnock, Testing Supervisor

End of Report

THIS REPORT SHALL NOT BE REPRODUCED EXCEPT IN ITS ENTIRETY AND WITH THE WRITTEN PERMISSION OF ELISA TECHNOLOGIES, INC.

2501 NW 66th Court, Gainesville, Florida 32653, USA
Telephone: (352) 337-3929 Fax: (352) 337-3928 E-mail: info@elisa-tek.com