

**ELISA Technologies, Inc.** 



elisa-tek.com

PAGE 1 OF 1 - REPORT 2205028E (of A-F)

**REPORT DATE:** May 13, 2022

REPORT TO: Eva Lau Mayway Corporation 1338 Mandela Parkway Oakland, CA 94607 Tel: 510-208-3113 x8131 Email: evalau@mayway.com

**SUBJECT:** PO Number: FU+JIU/CU/QU; Samples for Gluten Analysis

**RECEIVED:** May 4, 2022; UPS Ground TRK# 12 89E 637 03 6246 7251

## ANALYTICAL RESULTS: RESULTS APPLY ONLY TO THE SAMPLES AS RECEIVED.

CONTROL #	SAMPLE DESCRIPTION	ANALYSIS IDENTIFICATION	Gluten
2205028-40	MW# 3645 Fortify Stomach Teapills Xiang Sha Yang Wei Wan Batch# A9H56	500821 Quantitative <b>Gluten</b> ELISA	Not Detected

**NOTE:** The ALLER-TEK<sup>®</sup> Gluten assay is a sandwich ELISA for the quantitative analysis of gluten in raw materials and finished food products. Results are expressed in ppm gluten with a quantifiable range of 2.5-80 ppm; it detects the gliadin and glutenin fractions of gluten. A result of "not detected" indicates that gluten was not detected in the sample (<2.5 ppm). Recovery may be reduced in samples that have been hydrolyzed.

The ALLER-TEK  $^{\odot}$  Gluten ELISA has been validated and certified with license no. 081202 as a Performance Tested Method by the AOAC Research Institute.

Information in italics has been supplied by the customer.

## Respectfully Submitted for ELISA Technologies, Inc.

Ву

shanshand By

Shanshan Lin, R&D Supervisor

End of Report

THIS REPORT SHALL NOT BE REPRODUCED EXCEPT IN ITS ENTIRETY AND WITH THE WRITTEN PERMISSION OF ELISA TECHNOLOGIES, INC.

2501 NW 66<sup>th</sup> Court, Gainesville, Florida 32653, USA +1 352 337 3929 info@elisa-tek.com