



ELISA Technologies, Inc.

www.elisa-tek.com



PAGE 1 OF 1 - REPORT 1807158

REPORT DATE: July 20, 2018

REPORT TO: Eva Lau
Mayway Corporation
1338 Mandela Parkway
Oakland, CA 94607
Tel: 510-208-3113
Email: evalau@mayway.com

SUBJECT: Samples for Gluten Analysis; RUSH ORDER

RECEIVED: July 19, 2018; UPS Next Day Air TRK# 1Z 89E 637 01 4722 3135

ANALYTICAL RESULTS:

CONTROL #	SAMPLE DESCRIPTION	ANALYSIS IDENTIFICATION	Gluten
1807158-1	Item# 3368 Salvia Teapills/ Dan Shen Yin Wan Batch# 2016061	500821 Quantitative Gluten ELISA	30.3 ppm
██████████	██ ██	████████████████████ ████████████████████	████████████████ ████████████████

RESULTS APPLY ONLY TO THE SAMPLES TESTED

NOTES:

The ALLER-TEK Gluten assay utilizes ELISA techniques and is designed for the measurement of cereal gluten at low levels, nominally within a range of 5-80 PPM, in raw materials and finished food products; it detects the gliadin and glutenin fractions of gluten.

For the quantitation of barley gluten, a standard curve derived from barley proteins is used.

The ALLER-TEK Gluten ELISA has been validated and certified with license no. 081202 as a Performance Tested Method by the AOAC Research Institute.

Respectfully Submitted for **ELISA Technologies, Inc.**

By Tom d'angelo
Analyst

By Emily Warnock
Emily Warnock, Testing Supervisor

End of Report

THIS REPORT SHALL NOT BE REPRODUCED EXCEPT IN ITS ENTIRETY AND WITH THE WRITTEN PERMISSION OF ELISA TECHNOLOGIES, INC.

2501 NW 66th Court, Gainesville, Florida 32653, USA
Telephone: (352) 337-3929 Fax: (352) 337-3928 E-mail: info@elisa-tek.com