



ELISA Technologies, Inc.

www.elisa-tek.com



PAGE 1 OF 1 - REPORT 1807090

REPORT DATE: July 16, 2018

REPORT TO: Kelly Lau
Mayway Corporation
1338 Mandela Parkway
Oakland, CA 94607
Tel: 510-208-3113
Email: evalau@mayway.com

SUBJECT: PO Number: KL07062018; Samples for Barley Gluten Analysis

RECEIVED: July 11, 2018; UPS 3 Day Select TRK# 1Z 89E 637 12 5366 0989

ANALYTICAL RESULTS:

CONTROL #	SAMPLE DESCRIPTION	ANALYSIS IDENTIFICATION	Gluten
████████	████████████████████ ████████████████████	████████ ██████████ ████████████████████	████████
████████	████████████████████ ████████████████████ ████████████████████	████████ ██████████ ████████████████████	████████
1807090-3	Product Code # 3676 Xiao Yao Teapills Batch 201702100	500822 Quantitative Barley Gluten ELISA	> 80.0 ppm
████████	████████████████████ ████████████████████	████████ ██████████ ████████████████████	████████
████████	████████████████████ ████████████████████ ████████████████████	████████ ██████████ ████████████████████	████████

RESULTS APPLY ONLY TO THE SAMPLES TESTED

NOTES:

The ALLER-TEK Gluten assay utilizes ELISA techniques and is designed for the measurement of cereal gluten at low levels, nominally within a range of 5-80 PPM, in raw materials and finished food products; it detects the gliadin and glutenin fractions of gluten.

For the quantitation of barley gluten, a standard curve derived from barley proteins is used.

The ALLER-TEK Gluten ELISA has been validated and certified with license no. 081202 as a Performance Tested Method by the AOAC Research Institute.

Respectfully Submitted for **ELISA Technologies, Inc.**

By Emily [Signature]
Analyst

By Hayley A. Curilla [Signature]
Hayley Curilla, Production Supervisor

End of Report

THIS REPORT SHALL NOT BE REPRODUCED EXCEPT IN ITS ENTIRETY AND WITH THE WRITTEN PERMISSION OF ELISA TECHNOLOGIES, INC.

2501 NW 66th Court, Gainesville, Florida 32653, USA
Telephone: (352) 337-3929 Fax: (352) 337-3928 E-mail: info@elisa-tek.com